


Appendix 2


Example of Vest Pictorial




Criss-cross vests must always cross in front.

- **Left - Correct Application**
- **Right - Incorrect Application**

Examples of Bed Pictorials



Straps should always be snug, but should not interfere with breathing. You should be able to slide your open hand (flat) between the device and the patient. If too loose, the patient could climb over, under, around, through, or between the siderails and become suspended in the device, resulting in chest compression, strangulation and death. Do not criss-cross straps behind the back. The patient could turn around in bed, loosen the straps, and try to climb over the siderails.



Secure waist straps to the movable part of the bed frame adjacent to the patient's hips. NEVER attach the straps to the siderails. If there is a danger or a history of the patient sliding up or down the mattress, a product with a pelvic piece may be indicated - check with the proper medical authority. Sliding up or down may result in chest compression or suffocation.



Appendix 4

Examples of Wheelchair Pictorials

Products must be comfortable and not restrict breathing in any way. At the first sign of discomfort, consult with the proper medical authority to correct the problem or to find another solution. If the patient has or develops a history of sliding, a pelvic piece may be required.

Proper body alignment is important. The belt should be at a 45 degree angle to hold hips back in the chair. The belt should be snug but comfortable over the thighs to help minimize the possibility of the patient sliding. Adjust footrests so the feet are supported—improper support causes unnecessary pressure behind the knees and compromises blood circulation.

Do not use back cushions—if dislodged, they may allow excess slack in the straps. Use extra vigilance when patients use seat cushions—if dislodged, some patients may be able to slide down, resulting in serious injury or death.



Products that tie to prevent sliding should always be attached under the chair out of the patient's reach—NEVER tied behind the back

Examples of Sizing Guides*

ADULT SIZING GUIDE					
Binding Color	Size	Patient Weight		Chest Size	
		Lb	Kg	In.	Cm
White	XS	60-90	27-41	28-34	71-86
Red	S	80-120	36-54	32-39	81-99
Green	M	110-155	50-70	37-45	94-114
Yellow	L	145-190	66-86	43-49	109-124
Blue	XL	180-230	82-104	47-52	119-132
Black	XXL	220-275	100-125	50-56	127-142

TODDLER SIZING GUIDE		
Color	Size	Age
Red	Small	2 - 4 years
Green	Medium	4 - 6 years
Blue	Large	6 - 8 years

The charts provide general sizing guidelines only--the patient's individual physical characteristics may require a different size product. The proper medical authority must make the final decision on what size product to use.

FDA strongly recommends that the above referenced color codes for sizing be used in order to make color coding more universal.

* All pictorials and sizing guides adapted with permission from the Posey Co. Use of this information does not imply an endorsement of the Posey Co. by FDA.